

Instruction
For
Neitz Otoscope FO

1. Switch (See Drawing.)

- Sliding down the Slide SW ⑥ to turn on the light(The indication at the slide shows in red.)
- Slide up the same to turn off. (The indication is in white.)

2. Specula

Specula③ has two sizes of tips(2.5mm, 4mm). Choice the size to meet the patient's ears. When attaching the Specula to the Otoscope, rotate the Specula clock-wise with adjusting the projection④ in the Specula to the cut⑤ in the Otoscope. Please pay the attention that the projection④ could be damaged if the Specula is attached without the above adjusting.

3. Loupe

When observing objects with magnifying, use the loupe①. When the loupe is not necessary under medical treatment, turn it to be kept upwards.

4. Bulb Replacement (See Drawing.)

Remove the Otoscope head from Battery handle ⑦ with rotating it counter-clockwise. Light bulb gets a very hot, when it just burn out. Replace the light bulb after it was cooled enough. Remove the bulb from the Otoscope and insert new bulb ②. Use Neitz Spare light bulb No. L-69.

5. Interchange Batteries. (See Drawing.)

Remove the Otoscope head from Battery handle ⑦ with rotating it counter-clockwise. Take out the both of AA Cell Batteries and insert the new ones with its Minus (-) terminal fast. We recommend using Alkaline Cell Battery.

NEITZ INSTRUMENTS CO., LTD.

4F Ichibancho Court, 15-21, Ichibancho, Chiyoda-ku, Tokyo 102-0082, Japan

TEL: +81-3-3237-0552

