HUVITZ APPLANATION TONOMETER
HT-5000



Huvitz Re:define, Re+create



HUVITZ APPLANATION TONOMETER HT-5000

- *The world standard Goldmann applanation tonometer*
- *The digital display makes it easy to read even in a dark room.*
- *The wide grip comfortable knob for easier adjustment of measuring prism*
- *Optimized algorithm realizes the fast process of IOP calculation.*
- *Instantaneous wake-up time from sleep mode!*
- *The power supply by the most commonly used AAA batteries guarantees easy installation and maintenance.*
- *The slanted display matches to the angle of a user' s view for more comfortable observation.*

SPECIFICATIONS

Measuring force	generated by leverage weight
Installation	Mountable on peg on microscope.
Measuring range	3~75mm Hg
Reverse span	=/< 0.49 mN
Power Consumption	DC3V
Net weight	0.50kg (without accessories, with batteries)

Designs and details can be changed without prior notice for the purposes of improvement.

HUVITZ Co., Ltd. 38, Burim-ro 170beon-gil, Dongan-gu,
Anyang-si, Gyeonggi-do, 14055, Republic of Korea
Tel:+82-31-428-9100 Fax:+82-31-477-9022(F/A) http: //www.huvitz.com

(M1XXCL-22-00001, Rev A. 170426)